



Nitrogen
Determination
by Kjeldahl
Method

PanReac 
AppliChem
ITW Reagents

Determination of Nitrogen according to Kjeldahl



For longer than 100 years the Kjeldahl method has been used for the determination of nitrogen in a wide range of samples. The determination of Kjeldahl nitrogen is made in foods and drinks, meat, feeds, cereals and forages for the calculation of the protein content. Also the Kjeldahl method is used for the nitrogen determination in wastewaters, soils and other samples.

It is an official method and it is described in different normatives such as AOAC, USEPA, ISO, Pharmacopeias and different European Directives. The Kjeldahl method is used to determine the nitrogen content in organic and inorganic samples.

The procedure involves three major steps:

- **Digestion:** The sample is digested in boiling concentrated sulfuric acid, with the addition of a catalyst, until complete dissolution and oxidation. The nitrogen contained in the sample becomes Ammonium Sulfate.




- **Distillation:** Adding an excess of sodium hydroxide solution, the ammonium ion is released in ammonia form, distilled and received on a boric acid solution or a sulfuric or hydrochloric acid volumetric solution.



- **Titration:** The ammonia is determined with a volumetric acid solution or by back-titration with sodium hydroxide solution of a known concentration if it was received on hydrochloric or sulfuric acid. The results can be expressed in % N, % NH₃ or protein (%N x factor).

Kjeldahl equipment



The following table shows our Kjeldahl program. If you wish to receive more information, don't hesitate to request our General Catalogue or visit our web in www.itwreagents.com.

Code	Description	Package									
Digestion											
Catalyst											
173350	Kjeldahl Catalyst (Cu) (0.3% in CuSO ₄ ·5H ₂ O) tablets (Potassium Sulfate + Copper(II) Sulfate). Missouri Catalyst.		1.25 kg	(250 tablets of 5.0 g)							
			3.5 kg	(1000 tablets of 3.5 g)							
			5 kg	(1000 tablets of 5.0 g)							
174428	Kjeldahl Catalyst (Cu) (6.25% in CuSO ₄ ·5H ₂ O) tablets according to Directive 93/28/EEC (Potassium Sulfate + Copper(II) Sulfate)		4 kg	(1000 tablets of 4.0 g)							
175639	Kjeldahl Catalyst (Cu) (9% in CuSO ₄ ·5H ₂ O) tablets (Potassium Sulfate + Copper(II) Sulfate)		1650 g	(1000 tablets of 1.65 g)							
			5 kg	(1000 tablets of 5.0 g)							
172429	Kjeldahl Catalyst (Cu-Se) powder (Potassium Sulfate + Copper(II) Sulfate + Selenium). For N determination according to Wieninger.		1000 g	(1000 tablets of 1.0 g)							
			5 kg	(1000 tablets of 5.0 g)							
172926	Kjeldahl Catalyst (Cu-Se) (1.5% CuSO ₄ ·5H ₂ O + 2% Se) tablets Potassium Sulfate + Copper(II) Sulfate + Selenium). For N determination according to Wieninger.		1000 g	(1000 tablets of 1.0 g)							
			3.5 kg	(1000 tablets of 3.5 g)							
			5 kg	(1000 tablets of 5.0 g)							
175570	Kjeldahl Catalyst (Cu-Se) (9% CuSO ₄ ·5H ₂ O + 0.9% Se) tablets (Potassium Sulfate + Selenium metal + Copper(II) Sulfate 5-hydrate). For soil analysis.		4 kg	(1000 tablets of 4.0 g)							
173349	Kjeldahl Catalyst (Cu-TiO ₂) tablets (Potassium Sulfate + Sodium Sulfate + Copper(II) Sulfate 5-hydrate + Titanium(IV) Oxide)		875 g	(250 tablets of 3.5 g)							
			1.25 kg	(250 tablets of 5.0 g)							
			3.71 kg	(1000 tablets of 3.71 g)							
			5 kg	(1000 tablets of 5.0 g)							
173348	Kjeldahl Catalyst (Se) tablets (Potassium Sulfate + Selenium)		3.5 kg	(1000 tablets of 3.5 g)							
			5 kg	(1000 tablets of 5.0 g)							
177033	Kjeldahl Catalyst (Cu) (1.96% in CuSO ₄ ·5H ₂ O) tablets (Potassium Sulfate + Copper(II) Sulfate 5-hydrate)		5 kg	(1000 tablets of 5.0 g)							
177040	Kjeldahl Catalyst (Cu) (10.26% in CuSO ₄ ·5H ₂ O) tablets (Potassium Sulfate + Copper(II) Sulfate).		4 kg	(1000 tablets of 4.0 g)							
Acids and oxidants											
173163	Sulfuric Acid 98%		1000 mL		2.5 L		25 L				
121076	Hydrogen Peroxide 30% w/v (100 vol.) for analysis		1000 mL		5 L						
Antifoaming											
211628	Silicone antifoaming liquid		100 mL		250 mL		500 mL				
Distillation											
Alkalis for liberation of ammonia											
131687	Sodium Hydroxide pellets for analysis, ACS, ISO		500 g		1 kg		5 kg		25 kg		
141571	Sodium Hydroxide solution 50% w/v pure		5 L		10 L						
171220	Sodium Hydroxide solution 40% w/w		1000 mL		5 L		10 L		25 L		
122666	Sodium Hydroxide solution 32% w/v for analysis		1000 mL		2.5 L		5 L		10 L		25 L
Solutions for capture of ammonia											
282972	Boric Acid solution 1% for volumetric analysis		5 L								
287096	Boric Acid solution 2% for volumetric analysis		5 L		25 L						
282928	Boric Acid solution 3% for volumetric analysis		1000 mL								
282222	Boric Acid solution 4% for volumetric analysis		1000 mL		5 L		25 L				
283334	Ammonia Fixative solution 1% for volumetric analysis		5 L		25 L						
181023	Hydrochloric Acid 0.1 mol/l (0.1N) volumetric solution		1000 mL		2.5 L		5 L		10 L		
181022	Hydrochloric Acid 0.5 mol/l (0.5N) volumetric solution		1000 mL		5 L		10 L				
181061	Sulfuric Acid 0.05 mol/l (0.1N) volumetric solution		1000 mL		5 L		10 L				
181060	Sulfuric Acid 0.25 mol/l (0.5N) volumetric solution		1000 mL		2.5 L		10 L				
Titration											
Volumetric solutions											
181693	Sodium Hydroxide 0.1 mol/l (0.1N) volumetric solution		1000 mL		2.5 L		5 L		10 L		
181023	Hydrochloric Acid 0.1 mol/l (0.1N) volumetric solution		1000 mL		2.5 L		5 L		10 L		
181061	Sulfuric Acid 0.05 mol/l (0.1N) volumetric solution		1000 mL		5 L		10 L				
Indicators											
283303	Indicator 4.8, Mixed (Methyl Red-Bromocresol Green) for volumetric analysis		250 mL								
282430	Indicator 4.4, Mixed (Methyl Red-Methylene Blue) for volumetric analysis		250 mL								
281618	Methyl Red solution 0.1% for volumetric analysis		250 mL								

Packaging Symbols:



Polyethylene bottle



Polyethylene container (cubitainer) with tap



Polypropylene bucket with handle



Polyethylene canister



Glass bottle

Because of the continuous updating of our catalog of products and presentations, kindly confirm product availability, through our website.



A173,EN;201509

PanReac 
AppliChem
ITW Reagents

PanReac Química SLU

C/ Garraf 2, Polígono Pla de la Bruguera
E-08211 Castellar del Vallès
(Barcelona) Spain
Phone +34 937 489 400
Fax +34 937 489 401
central@panreac.com

AppliChem GmbH

Ottoweg 4
DE-64291 Darmstadt
Germany
Phone +49 6151 9357 0
Fax +49 6151 9357 11
service@de.applichem.com

Nova Chimica Srl

Via G. Galilei, 47
I-20092 Cinisello Balsamo
(Milano) Italy
Phone +39 02 66045392
Fax +39 02 66045394
info@novachimica.com

www.itwreagents.com