

Specification

AppliFect LowTox

A9027

Synonym	Transfection reagent for very sensitive cells
additional product description	Reagent for transfection of very sensitive, eukaryotic cells.
HS-No.:	38220000
Storage:	-20°C
WGK:	<ul style="list-style-type: none"> • lowest toxicity • very high expression rate • approx. 500 (24-well) or approx. 200 (6-well) transfections per 1 ml nwg
<p>Comment</p> <p>AppliFect LowTox is an improvement of AppliFect. While developing this new reagent, the goal was to achieve highest compatibility with many cell lines and highest transfections efficiency. The transfected DNA is rapidly released within the cytosol from the lipoplex of cationic lipids and colipids. The reagent is supplied as a sterile, aqueous solution. AppliFect LowTox is shipped at ambient temperature. Please freeze upon arrival.</p>	