


## Specification

**AppliFect**

**A8886**

<b>Synonym</b>	Transfection reagent AppliFect
<b>additional product description</b>	Reagent for transfection of eukaryotic cells with nucleic acids
<b>HS-No.:</b>	38220000
<b>Storage:</b>	2-8°C
	<ul style="list-style-type: none"> <li>• for transient or stable transfection</li> <li>• for adherent or suspension cells</li> <li>• no serum inhibition</li> <li>• applicable for sensitive cells</li> <li>• polycationic reagent</li> <li>• sufficient for approx. 200 (6-well) transfections per ml</li> </ul>
<b>WGK:</b>	nwg
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • [info.de@itwreagents.com](mailto:info.de@itwreagents.com) • [www.itwreagents.com](http://www.itwreagents.com)

CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900

IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340