

Specification

PBS buffer (1X, Dulbecco's) - Powder**A0964**

Synonym	Powder mixture fo prepare phosphate buffered saline
	9.55
HS-No.:	38210000
Storage:	RT
LGK:	10 - 13
WGK:	1
	hygroscopic
Specification	
Osmolality	262 - 290 mOsm/kg H ₂ O
pH (25°C)	7.2 - 7.8
Composition:	
di-Sodium hydrogen phosphate anhydrous	1150.00 mg/L
Potassium chloride	200.00 mg/L
Potassium dihydrogen phosphate	200.00 mg/L
Sodium chloride	8000.00 mg/L